



People.Health.Care.

Tribute® Wrap Full Leg **NEW**

One-piece coverage. One-of-a-kind adjustability.

Proven foam technology

Promotes a micro-massage-type effect to break up and soften hardened skin or fibrotic tissue to support skin health

Chevron channeling

Designed based on MLD principles of directing fluid to collateral pathways

Easy-to-pull finger grip

Conforms for the right, comfortable fit to support self-management

Supports mobility

Angled straps at flexion points for a comfortable and functional design



More options, more coverage

40 unique combinations of size, length, orientation and width to accommodate more patients.



Cost-conscious

Reimbursable under HCPCS code A6526 to improve patient access.





People.Health.Care.

Tribute® Wrap Full Leg

Off-the-shelf adjustable foam gradient compression garment

Tribute® Wrap Full Leg supports easy self-care of lymphedema and edema of the leg and foot during the evening, night, or low activity.

Tribute Wrap comfortably adjusts to your unique shape to help maintain gains made during therapy, while allowing the flexibility to support further reduction.

Indications

- Various stages of lymphedema
- Edema



Black exterior
Slate interior

- 6-month warranty
- Easy garment care in washer and dryer

Tribute® Wrap Full Leg

left or right orientation, sold individually, Black Sleep Sleeve included

Size	S	S WIDE	M	M WIDE	L	L WIDE	XL	XL WIDE	XXL	XXL WIDE	
G°	53–60 cm	59–66 cm	59–66 cm	65–72 cm	65–72 cm	71–78 cm	71–78 cm	77–84 cm	77–84 cm	83–90 cm	
E°	42–48 cm	47–53 cm	47–53 cm	52–58 cm	52–58 cm	57–63 cm	57–63 cm	62–68 cm	62–68 cm	67–73 cm	
D°	31–36 cm		35–40 cm		39–44 cm		43–48 cm		47–53 cm		
C°	30–36 cm		33–41.5 cm		36–45 cm		40.5–49 cm		44–54 cm		
B°	24–30 cm		27–33 cm		30–36 cm		34–40 cm		37–43 cm		

Length	regular	long
AY ^L	up to 22 cm	up to 26 cm
YG ^L	up to 80 cm	up to 88 cm

How to Order

Tribute® Wrap Full Leg

1 Measure

Only 7 measurements needed.

2 Order

Available via your preferred garment supplier:

- Local DME
- National provider
- Online distributor

[US](http://www.lohmann-rauscher.com/us-en/where-to-buy) www.lohmann-rauscher.com/us-en/where-to-buy

[CA](http://www.lohmann-rauscher.ca) www.lohmann-rauscher.ca

Dealers can purchase direct from L&R via phone, fax, or online store.